

# Pro 30 Laboratory Reciprocal Shaker

## Instruction Manual

**Catalog Numbers:**

S2030-RC-B

S2030-RC-B-230V



## Table of Contents

1.0 Introduction .....	1
2.0 Packing .....	1
3.0 Specifications .....	1
4.0 Safety Information .....	1
5.0 Installation .....	2
5.1. Device Placement .....	2
5.2. Attaching the Power Cord .....	2
5.3. Mounting Platform .....	2
6.0 Operating Instructions .....	2
6.1. Control Panel .....	2
7.0 Troubleshooting .....	3
8.0 Maintenance and Cleaning .....	3
8.1. Cleaning .....	3
9.0 Installation of Double Platform .....	4
10.0 Limited Warranty .....	5
11.0 Equipment Disposal .....	5

## 1.0 Introduction

Congratulations on purchasing your Labnet Pro 30 Laboratory Reciprocal Shaker. This shaker features a left to right linear action to provide gentle to vigorous shaking. Four platform configurations make this model valuable for a variety of applications.

## 2.0 Packing

The Labnet Pro 30 Laboratory Reciprocal Shaker includes the following components:

- ▶ Pro 30 Reciprocal Shaker
- ▶ Instruction Manual

**NOTE:** Platform needs to be purchased separately.

## 3.0 Specifications

<b>Electric Power Supply</b>	115V ± 10%, 50/60Hz, 0.2A or 230V ± 10%, 50/60Hz, 0.1A 230V includes EU and UK power cords
<b>Fuse</b>	2 x 0.5A (115V); 2 x 0.25A (230V)
<b>Speed</b>	From 5 to 150 strokes per minute
<b>Stroke Length</b>	30 mm
<b>Timer</b>	20 to 120 minutes, in 10 minute steps, HOLD function
<b>Load</b>	11 lbs. (5 kg) including platform
<b>Dimensions (W x D x H)</b>	10 x 12.6 x 6.3 in. (25.5 x 32 x 16 cm)
<b>Platform Dimensions (W x D)</b>	12 x 12 in. (30 x 30 cm)
<b>Weight</b>	15 lbs. (6.6 kg)
<b>Operating Conditions</b>	4°C to 70°C, up to 85% RH, non-condensing

## 4.0 Safety Information

- ▶ Do not use the shaker near water sources. Be sure that water doesn't get on the device, especially during the cleaning procedures.
- ▶ Do not expose the shaker to aggressive vapors (e.g. organic solvents).
- ▶ Do not insert objects into the shaker (such as screw driver or metal rods).
- ▶ Never operate the shaker without platform securely attached.
- ▶ Hazardous samples should be shaken only in appropriate containment vessels.

**NOTE:** If the equipment is used in a manner not specified by the manufacturer, the protection provided by the equipment may be impaired.



**Do not shake flammable or explosive samples.**

## 5.0 Installation

### 5.1. Device Placement

When selecting the right place for device, please consider following:

- ▶ Put the device on smooth, horizontal and stable place.
- ▶ Leave enough space beyond the device for normal air circulation, min. 10 cm.
- ▶ Leave enough space around the device, for easy access and maintenance.
- ▶ Avoid locations subject to shocks and vibrations.
- ▶ Avoid locations subject to frequent power fluctuations or power losses, it may affect shaker performance.

**NOTE:** There has to be easy access to the shaker control panel and main plug in case of emergency.

### 5.2. Attaching the Power Cord

First connect the main power cord to the shaker then connect the plug end of the cord to a grounded wall socket. Avoid lines on which powerful electric motors, refrigerators and similar devices are connected.



- ▶ Be careful when you plug the cord to a grounded wall socket.
- ▶ Do not touch the plug with wet hands.
- ▶ Do not pull the plug by the cord.

### 5.3. Mounting Platform

Mount the accessory platform by firmly inserting the four platform mounting pins into the four rubber grommet points on the shaker.

## 6.0 Operating Instructions

### 6.1. Control Panel



1. **POWER switch** – to turn the shaker ON (power switch illuminates) or OFF.
2. **SPEED potentiometer** – to set the shaker speed.
3. **TIME potentiometer** – to set the timed operation.
4. **HOLD potentiometer** – yellow light illuminates when the shaker is set to HOLD (continuous operation).
5. **START/STOP button** – to START or STOP the shaking operation.
6. **RUN** – green light illuminates when the shaker is running.

**NOTE:** The speed can be adjusted at any time during shaker operation. The time, however, can not be changed while the shaker is in motion. To change the time press the STOP button first, allow the shaker platform to come to a full stop and set the time value again. Restart the shaker.

## 7.0 Troubleshooting

Should you have a question about the operation of the Labnet Pro 30 Reciprocal Shaker or if service is required, contact your Corning representative. Do not send in a unit for service without first calling to obtain a repair authorization number. Should the unit require return to Corning for service, it should be properly packed to avoid damage. Any damage resulting from improper packaging shall be the responsibility of the user.

Problem	Explanation/Solution
POWER switch doesn't illuminate when turned ON.	<ul style="list-style-type: none"><li>▶ Check the connection to main power supply.</li><li>▶ Check fuses</li><li>▶ Call authorized service</li></ul>
Shaker doesn't start when you press START/STOP button.	<ul style="list-style-type: none"><li>▶ Check the connection to main power supply.</li><li>▶ Check fuses</li><li>▶ Call authorized service</li></ul>
HOLD yellow signal light and RUN green signal light are blinking.	<ul style="list-style-type: none"><li>▶ Motor overload - reduce the load on the platform.</li><li>▶ Call authorized service.</li></ul>

## 8.0 Maintenance and Cleaning

- ▶ No routine maintenance or lubrication is required.
- ▶ It is recommended that the shaker is kept clean to ensure long trouble-free operation.

### 8.1. Cleaning

The housing of the shaker can be cleaned with a soft, damp cloth or special cleaners for plastic surfaces (e.g. for car interior – armature). Do not immerse the device or spill liquids on the device.


**NOTE:** Do not use any aggressive or abrasive cleaners (acetone, nitro, polish etc.) because the surface can get permanently damaged.

 Prior to any cleaning or maintenance of the device, unplug the power cord from the wall socket.

## 9.0 Installation of Double Platform

Before the installation of double platform, turn off the POWER switch.

For installation of double platform, prepare upper platform with rubber pad and posts kit, which includes four posts, four washers, four screws and four nuts.



1. Set four posts (A) above four holes at the corners of the bottom platform.
2. Put the washer (B) on each screw (C) and then screw it into each post from the bottom side of the platform.
3. Fit upper platform (D) on the top of four posts, so that the threads of the posts are protruding through four corner holes.
4. Screw four nuts (E) on the threads of four posts, to fix the upper platform.

## 10.0 Limited Warranty

Corning Incorporated (Corning) warrants that this product will be free from defects in material and workmanship for a period of one (2) years from date of purchase. CORNING DISCLAIMS ALL OTHER WARRANTIES WHETHER EXPRESSED OR IMPLIED, INCLUDING ANY IMPLIED WARRANTIES OF MERCHANTABILITY OR OF FITNESS FOR A PARTICULAR PURPOSE. Corning's sole obligation shall be to repair or replace, at its option, any product or part thereof that proves defective in material or workmanship within the warranty period, provided the purchaser notifies Corning of any such defect. Corning is not liable for any incidental or consequential damages, commercial loss or any other damages from the use of this product.

This warranty is valid only if the product is used for its intended purpose and within the guidelines specified in the supplied instruction manual. This warranty does not cover damage caused by accident, neglect, misuse, improper service, natural forces or other causes not arising from defects in original material or workmanship. This warranty does not cover motor brushes, fuses, light bulbs, batteries or damage to paint or finish. Claims for transit damage should be filed with the transportation carrier.

In the event this product fails within the specified period of time because of a defect in material or workmanship, contact Corning Customer Service at: USA/Canada 1.800.492.1110, outside the U.S. +1.978.442.2200, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences), or contact your local support office.

Corning's Customer Service team will help arrange local service where available or coordinate a return authorization number and shipping instructions. Products received without proper authorization will be returned. All items returned for service should be sent postage prepaid in the original packaging or other suitable carton, padded to avoid damage. Corning will not be responsible for damage incurred by improper packaging. Corning may elect for onsite service for larger equipment.

Some states do not allow limitation on the length of implied warranties or the exclusion or limitation of incidental or consequential damages. This warranty gives you specific legal rights. You may have other rights which vary from state to state.

No individual may accept for, or on behalf of Corning, any other obligation of liability, or extend the period of this warranty.

For your reference, make a note of the serial and model number, date of purchase, and supplier here.

Serial No. \_\_\_\_\_ Date Purchased \_\_\_\_\_  
Model No. \_\_\_\_\_ Supplier \_\_\_\_\_

## 11.0 Equipment Disposal



According to Directive 2012/19/EU of the European Parliament and of the Council of 4 July 2012 on waste electrical and electronic equipment (WEEE), this product is marked with the crossed-out wheeled bin and must not be disposed of with domestic waste.

Consequently, the buyer shall follow the instructions for reuse and recycling of waste electronic and electrical equipment (WEEE) provided with the products and available at [www.corning.com/weee](http://www.corning.com/weee).

To request certificates, please contact us at [www.labnetlink.com](http://www.labnetlink.com).

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For additional product or technical information, visit [www.corning.com/lifesciences](http://www.corning.com/lifesciences) or call 800.492.1110. Outside the United States, call +1.978.442.2200 or contact your local Corning sales office.

**CORNING**

Corning Incorporated  
Life Sciences

[www.corning.com/lifesciences](http://www.corning.com/lifesciences)

**NORTH AMERICA**

t 800.492.1110  
t 978.442.2200

**ASIA/PACIFIC**

**Australia/New Zealand**  
t 61 427286832

**Chinese Mainland**  
t 86 21 3338 4338

**India**

t 91 124 4604000

**Japan**

t 81 3-3586 1996

**Korea**

t 82 2-796-9500

**Singapore**

t 65 6572-9740

**Taiwan**

t 886 2-2716-0338

**EUROPE**

HTL.sales@corning.com  
HTL.service@corning.com

**LATIN AMERICA**

grupoLA@corning.com

**Brazil**

t 55 (11) 3089-7400

**Mexico**

t (52-81) 8158-8400