

Labnet™ BioPette™ Plus Single- and Multi-channel Pipettors

Engineered to provide a high level of accuracy, reproducibility, and comfort

Available in single-, 8- and 12-channel configurations, featuring smooth plunger movement with low pipetting forces. Volume can be set by turning either the thumb wheel or the plunger button and is protected by a locking mechanism. BioPette Plus pipettors have been designed to maintain their calibration, even after repeated use and autoclaving.

A calibration key is included with each pipettor for easy in-house calibration.

BioPette Plus multi-channel pipettors feature individually spring-loaded piston assemblies for positive tip sealing and a curved ejector bar requiring minimal ejection force. The multi-channel manifold rotates 360° for ejecting and pipetting in any position.

All BioPette Plus pipettors feature an accurate and precise pipetting mechanism, enclosed in an ergonomic, well balanced and robust handle.

Features

- Factory calibrated to strict ISO 8655 standards
- 3-year warranty
- Fully autoclavable and UV-resistant
- Designed for universal fit with all common brands of pipette tips
- Adjustable tip ejector



Ordering Information



Labnet™ BioPette™ Plus Single- and Multi-channel Pipettors

Cat. No.	Volume Range (µL)	Accuracy (%)	Precision (%)	Non-filtered Tip (µL)
P3942-2	0.1-2	±40.0 to ±1.5	≤12.0 to ≤0.7	10
P3942-10	0.5-10	±4.0 to ±0.5	≤2.8 to ≤0.4	10
P3942-20	2-20	±3.0 to ±0.8	≤1.5 to ≤0.3	200
P3942-100	10-100	±3.0 to ±0.8	≤1.0 to ≤0.2	200
P3942-200	20-200	±2.0 to ±0.6	≤0.7 to ≤0.2	200
P3942-1000	100-1,000	±2.0 to ±0.6	≤0.6 to ≤0.2	1,000
P3942-5000	500-5,000	±1.2 to ±0.5	≤0.5 to ≤0.15	5,000
P3942-10000	1,000-10,000	±2.5 to ±0.5	≤0.6 to ≤0.2	10,000
P4808-10	0.5-10	±10.0 to ±2.0	≤8.0 to ≤1.2	10
P4808-50	5-50	±4.0 to ±1.6	≤2.5 to ≤0.6	200
P4808-200	20-200	±3.0 to ±1.0	≤1.5 to ≤0.6	200
P4808-300	50-300	±1.6 to ±1.0	≤1.5 to ≤0.6	300
P4812-10	0.5-10	±10.0 to ±2.0	≤8.0 to ≤1.2	10
P4812-50	5-50	±4.0 to ±1.6	≤2.5 to ≤0.6	200
P4812-200	20-200	±3.0 to ±1.0	≤1.5 to ≤0.6	200
P4812-300	50-300	±1.6 to ±1.0	≤1.5 to ≤0.6	300

Labnet BioPette Plus Starter Kit

Cat. No.	Description
P3942-SK4	Labnet BioPette Plus Starter Kit Includes 4 single-channel pipettors with calibration keys: 0.5 to 10 µL (Cat. No. P3942-10), 2 to 20 µL (Cat. No. P3942-20), 20 to 200 µL (Cat. No. P3942-200), 100 to 1,000 µL (Cat. No. P3942-1000); universal carousel stand for six single-channel pipettors (Cat. No. P3985-P); 3 racks of pipette tips: 10 µL (Cat. No. 301-03-102), 200 µL (Cat. No. 301-02-152), and 1,000 µL (Cat. No. 301-01-302)

Accessories

Cat. No.	Description	Qty/Pk
SP19475	Shelf clip for one single- or multi-channel pipettor	1
P3985-P	Universal carousel stand for six single-channel pipettors (up to 1,000 µL)	1
SP19476	Shaft filter for 5,000 and 10,000 µL pipettors	10
SP19477	Shaft filter for 5,000 and 10,000 µL pipettors	50

Related Products

Cat. No.	Description	Qty/Pk
R8300	Storage box with hinged lid, 100 x 1.5/2.0 mL, natural	1
R8300-A	Storage box with hinged lid, 100 x 1.5/2.0 mL, assorted color pack	5
R8300-BLK	Storage box with hinged lid, 100 x 1.5/2.0 mL, black opaque	1



www.labnetlink.com