

Labnet FastPette™ V2 Pipet Controller

Lightweight and ergonomic design makes the V2 a perfect fit for any lab

The FastPette V2 is a battery powered pipet controller ideal for liquid handling with 0.5 to 100 mL glass or plastic serological pipets. Designed with the user in mind, the FastPette has a lightweight handle with smooth finger triggers and ergonomically located switches.

The aspiration speed can be set to high or low to adjust for the volume of the pipet and viscosity of the liquid. The suction can also be adjusted by the amount of pressure applied to the trigger button. Liquids can be dispensed by gravity or by blow out using the motorized pump.

The FastPette V2 is protected from the entrance of liquids by a PTFE filter. To protect samples from cross contamination, both the filter and pipette holder are easily exchanged and autoclaved.

Features

- ▶ Works with pipets from 0.5 mL to 100 mL
- ▶ Cordless work for up to 8 continuous hours
- ▶ Ergonomically designed
- ▶ Autoclavable pipet holder
- ▶ 1 year warranty





Ordering Information

Labnet FastPette™ V2 Pipet Controller

Cat. No.	Description	Qty/Cs
P2000	FastPette V2 is supplied with a universal power supply*, wall hanger, bench stand, and 0.2 µm PTFE filter	1

*Universal power supply is compatible for US, EU, UK and AU.

Accessories

Cat. No.	Description	Qty/Cs
SP9143	PTFE autoclavable filter, 0.2 µm	5
SP9144	PTFE autoclavable filter, 0.45 µm	5
SP9029	Wall hanger	1
SP19030	Bench stand	1
SP9022	Nosepiece	1
SP29054	Pipet holder	1