

Pipettor Selection Guide

Labnet offers a wide selection of manual and electronic pipetting products to choose from. Select your liquid handling product based on the features listed below.

Product	Volume Settings	Channels	Volume Range	Pipet Types	Electronic Operation	Finger hook	Auto-clavable*	Multi-Dispensing	In-Lab Calibration
BioPette™ Plus Single-channel Pipettors	Push button Adjustment knob	Single	0.1-10,000 µL	N/A	X	✓	✓	X	✓
BioPette Plus Multi-channel Pipettors	Push button Adjustment knob	8 or 12	0.5-300 µL	N/A	X	✓	✓	X	✓
BioPette A Single-channel Pipettors	Push button Adjustment knob	Single	0.1-10,000 µL	N/A	X	X	✓	X	✓
BioPette A Multi-channel Pipettors	Push button Adjustment knob	8 or 12	1-300 µL	N/A	X	X	✓	X	✓
Labpette™ FX Fixed Volume Pipettors	Fixed volume	Single	2-1,000 µL	N/A	X	✓	✓	X	✓
Excel™ Single-channel Electronic Pipettors	Programmable	Single	0.1-1.0 µL	N/A	✓	✓	✓	✓	X
Excel Multi-channel Electronic Pipettors	Programmable	8 or 12	0.5-1,200 µL	N/A	✓	✓	✓	✓	X
Labpette R Repeating Pipettors	Click stop	Single	1-5,000 µL	N/A	X	X	X	✓	N/A
FastPette™ V2 Pipet Controllers	Finger switches	Single	0.5-100 mL	Glass or plastic	✓	N/A	✓	✓	N/A
FastPette Pro Pipet Controllers	Finger switches	Single	0.5-100 mL	Glass or plastic	✓	N/A	✓	✓	N/A

*NOTE: Please read the instruction manual before device autoclaving.

