













Spectrafuge 6C Tube Compatibility Chart

| | | Standard Blood Collection Tubes | | | Sarstedt** Blood Collection Tubes | | | | | | | |
|---------------------|---|---|---|---|---|---|---|---|---|---|---|---|
| Dimensions (mm) | 17 x 120 | 16 x 100 | 13 x 100 | 12 x 75 | 15 x 102 | 16 x 92 | 15 x 92 | 13 x 90 | 11 x 92 | 11 x 66 | 13 x 65 | 8 x 66 |
| Capacity* | 15 mL | 10 mL | 7 mL | 5 mL | 10 mL | 9-10 mL | 7-8 mL | 4.9 mL | 4-5 mL | 2-3 mL | 2-3 mL | 1-1.4 mL |
| Tube image |  |  |  |  |  |  |  |  |  |  |  |  |
| Radius (cm) | 8.2 | 8.5 | 7.7 | 7.7 | 8.5 | 8.5 | 8.5 | 7.7 | 7.7 | 7.7 | 7.7 | 7.7 |
| Max. RFC | 3,873 | 4,015 | 3,637 | 3,637 | 4,015 | 4,015 | 4,015 | 3,637 | 3,637 | 3,637 | 3,637 | 3,637 |
| Tubes per Rotor | 6 | 6 | 6 | 6 | 6 | 6 | 6 | 6 | 6 | 6 | 6 | 6 |
| Adapter (C0200-17A) | — | — | Yes | Yes | — | — | — | Yes | Yes | Yes | Yes | Yes |

*Capacity may vary based on serum volume

**Sarstedt is a registered trademark of Sarstedt Company. The compatibility of these tubes has been independently determined by Corning and there should be no implication of sponsorship or endorsement by Sarstedt Company.